


HA (A/Japanese white-eye/Hong Kong/1038/2006 (H5N1))



Product Name	HA (A/Japanese white-eye/Hong Kong/1038/2006 (H5N1))
Catalog No	IT-003-00514p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA (A/Japanese white-eye/Hong Kong/1038/2006 (H5N1)) (Genebank No. ABJ96775) protein (a.a. 17-530) with deletion at its cleavage site (Δ RRRKKR). Majority of the protein exists as trimer/oligomer
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	$\geq 95\%$
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified HA (A/Japanese white-eye/Hong Kong/1038/2006 (H5N1)) protein