




HA(Δ TM)(H₃N₂/ Wisconsin/67/X-161/2005)

Product Name	HA(Δ TM)(H ₃ N ₂ / Wisconsin/67/X-161/2005)
Catalog No	IT-003-0041 Δ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM) (H ₃ N ₂ / Wisconsin/67/X-161/2005)(a.a.18-530)(Genebank No. ABO37609)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H₃N₂/ Wisconsin/67/X-161/2005)

Note: