




HA(Δ TM)(H₇N₇)

Product Name	HA(Δ TM)(H ₇ N ₇)
Catalog No	IT-003-008 Δ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(A/chicken/Netherlands/1/03(H7N7)) protein (a.a.17-527) (Genebank No. AY338459)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H₇N₇)

Note: