




Pp150 (UL32)(Cytomegalovirus)

Product Name	Pp150 (UL32)(Cytomegalovirus)
Catalog No	IT-005-022Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in 500mM NaCl, 50mM Phosphate buffer (pH 7.4), 200mM imidazole, 8M urea. Reconstitute the protein using 100 µl Ultrapure water.
Protein	C-terminal 6xHis tagged Pp150 (UL32) immunodominant region (a.a. 862-1048)
Applications	WB standard, antibody ELISA, etc
Purity	>90%
Storage	Store at -80°C; <u>AVOID FREEZING/THAWING</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Pp150 (UL32) protein (a.a. 862-1048)

Note: