



HA(Δ TM)(A/Ohio/01/2007(H₁N₁))

Product Name	HA(Δ TM)(A/Ohio/01/2007(H ₁ N ₁))
Catalog No	IT-003-SW2 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(A/Ohio/01/2007(H ₁ N ₁)) protein (a.a. 18-529) (Genebank No. ACQ73203). This protein is from a virus closely related to the current circulating “novel human H ₁ N ₁ ” influenza virus as described in “V. Shinde et al: Triple-Reassortant Swine Influenza A (H1) in Humans in the United States, 2005–2009 , published at www.nejm.org May 7, 2009 (10.1056/NEJMoa0903812).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

