


HA(Δ TM)(A/South Carolina/1/18/ (H₁N₁))



Product Name	HA(Δ TM)(A/South Carolina/1/18 (H ₁ N ₁))
Catalog No	IT-003-0014 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(A/South Carolina/1/18 (H ₁ N ₁)) protein (a.a. 19-529)(Genebank No. AAD17229).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(A/South Carolina/1/18 (H₁N₁))