

HA(Δ TM)(A/Japan/305/1957(H₂N₂))



Product Name	HA(Δ TM)(A/Japan/305/1957(H ₂ N ₂))
Catalog No	IT-003-0032 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (0.5 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(A/Japan/305/1957(H ₂ N ₂)) protein (a.a. 19-529)(Genebank No. AAA43185).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

