

HA(Δ TM)(A/chicken/Vietnam/NCVD-016/08 (H5N1))



Product Name	HA(Δ TM)(A/chicken/Vietnam/NCVD-016/08(H5N1))
Catalog No	IT-003-00518p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(A/chicken/Vietnam/NCVD-016/08(H5N1)) protein (a.a.18-534) (Genebank No. ACO07033)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(A/chicken/Vietnam/NCVD-016/08(H5N1))