

HA₁(A/Albany/12/1951(H₁N₁))



Product Name	HA ₁ (A/Albany/12/1951(H ₁ N ₁))
Catalog No	IT-003-0017p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA ₁ (A/Albany/12/1951(H ₁ N ₁)) protein (a.a.18-345) (Genebank #: ABP49481)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

