

gp120(HIV-1/Clade CRF02_AG/ 01CM_1475MV)



Product Name	gp120(HIV-1/Clade CRF02_AG/ 01CM_1475MV)
Catalog No	IT-001-155p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp120(HIV-1/Clade CRF02_AG/ 01CM_1475MV) protein (a.a.30-503)(Genebank #: AAW72255)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(HIV-1/Clade CRF02_AG/
01CM_1475MV)