


gH(DI-II)/gL/gp42 (Ectodomain)(EBV) (Strain B95-8) Complex



Product Name	gH(DI-II)/gL/gp42(Ectodomain)(EBV)(Strain B95-8) Complex
Catalog No	IT-005-033p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged ectodomain of gH (EBV) protein (a.a. 1-344) (Genebank #: P03231) in complex with gp42 ectodomain (EBV)(Strain B95-8) (Genebank #: YP_401672)(a.a. 31-223), and gL(a.a. 25-137) (Genebank #: CAD53428).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gH(DI-II)/gL/gp42(Ectodomain)(EBV)(Strain B95-8) Complex