


HA(Δ TM)(H₅N₁)(A/Hubei/1/2010)



Product Name	HA(Δ TM)(H ₅ N ₁)(A/Hubei/1/2010)
Catalog No	IT-003-00519 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged ectodomain of HA(H ₅ N ₁)(A/Hubei/1/2010) protein(a.a. 18-530)(Genebank #: EPI267032)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H₅N₁)(A/Hubei/1/2010)