

# HA<sub>1</sub>(H<sub>5</sub>N<sub>1</sub>)( A/Hubei/1/2010)



---

Product Name	HA <sub>1</sub> (H <sub>5</sub> N <sub>1</sub> )( A/Hubei/1/2010)
Catalog No	IT-003-00519p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA <sub>1</sub> (H <sub>5</sub> N <sub>1</sub> )( A/Hubei/1/2010) protein(a.a. 18-341)(Genebank #: EPI267032)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA<sub>1</sub>(H<sub>5</sub>N<sub>1</sub>)( A/Hubei/1/2010)