




HA (B/Malaysia/2506/2004)

Product Name	HA (B/Malaysia/2506/2004)
Catalog No	IT-003-0054p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.
Protein	6xHis tagged HA (B/Malaysia/2506/2004) protein (a.a. 12-547)
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified HA (B/Malaysia/2506/2004) protein