



Anti-gp120(Clade A, B, C)

Product Name	Anti-gp120(Clade A, B, C)
Catalog No	IT-001-002M10
Source	Mouse monoclonal antibody, clone 21D6
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant gp140(B) protein (a.a. 24-680)(IT-001-0021p) (Genebank accession M63929).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with certain HIV-1 Clade A, B, C gp120 proteins as tested so far. Minimal cross-reactivity to other subtype HIV-1 envelop proteins.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-001-002M10

