

Anti-gp120(HXBc2), (Du151)(C)



Product Name	Anti-gp120(HXBc2) , (Du151)(C)
Catalog No	IT-001-002M12
Source	Mouse monoclonal antibody, clone 78D7
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant gp140(B) protein (a.a. 24-680)(IT-001-0021p) (Genebank accession M63929).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with HIV-1 gp120(HXBc2) and gp120(Du151)(C) proteins as tested so far . Minimal cross-reactivity to other subtype HIV-1 envelop proteins.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-001-002M12

