

# Anti- Gag(SIV/mac239)



---

Product Name	Anti- Gag(SIV/mac239)
Catalog No	IT-001-024
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Recombinant Gag(SIV/mac239) protein (Genebank accession # M33262).
Applications	ELISA(1:1000-1:100,000), IP, etc.; NOT for WB.
Specificity	Reacts with SIV Gag protein.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.