


Anti-Urease subunit alpha [Helicobacter pylori J99]



Product Name	Anti-Urease subunit alpha [Helicobacter pylori J99]
Catalog No	IT-021-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Urease subunit alpha [Helicobacter pylori J99] (Genebank accession # NP_222789)
Applications	Western Blotting (1:200-1:2000), ELISA, etc.
Specificity	Reacts with Urease subunit alpha [Helicobacter pylori J99]. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Urease subunit alpha [Helicobacter pylori J99]

Note: