

# E2 (Sindbis virus)



---

Product Name	E2 (Sindbis virus)
Catalog No	IT-022-002Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	C-terminal 6xHis tagged Sindbis virus E2 protein (a.a. 137-493) (Genebank Accession No. ACU25469).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Stable for 6 months at -20°C from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

