

# NP (Rift Valley Fever Virus)



---

Product Name	NP (Rift Valley Fever Virus)
Catalog No	IT-017-002p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged NP (RVFV) protein (a.a. 1-245) (Genebank Accession No: P21700).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

