

Nucleocapsid protein N (Canine distemper virus)



Product Name	Nucleocapsid protein N(CDV)
Catalog No	IT-007-001Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water.
Protein	C-terminal 6xHis tagged nucleocapsid protein N [Canine distemper virus strain Onderstepoort] protein (a.a. 1-523) (Genebank Accession No. AAK54665).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

