


# NS1 (St. Louis encephalitis virus)



---

Product Name	NS1 (St. Louis encephalitis virus)
Catalog No	IT-006-0058p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal V5-6xHis tagged NS1 (St. Louis encephalitis virus) antigen (a.a.791-1145) (Genbank Accession No. ACB58159).
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>90% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NS1 (St. Louis encephalitis virus)