

# gH(DI-III)/gL/gp42 (Ectodomain)(EBV) (Strain B95-8) Complex



Product Name	gH(DI-III)/gL/gp42(Ectodomain)(EBV)(Strain B95-8) Complex
Catalog No	IT-005-032p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged ectodomain of gH (EBV) protein (a.a. 1-529) (Genebank #: P03231) in complex with gp42 ectodomain (EBV)(Strain B95-8) (Genebank #: YP_401672)(a.a. 31-223), and gL(a.a. 25-137) (Genebank #: CAD53428).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

E  
250 kD  
150 kD  
100 kD  
75 kD  
50 kD  
37 kD  
25 kD  
20 kD  
15 kD  
10 kD

E: gH(DI-III)/gL/gp42(Ectodomain)(EBV)(Strain B95-8) Complex

