



Pp28(UL99) (Cytomegalovirus)

Product Name	Pp28(UL99) (Cytomegalovirus)
Catalog No	IT-005-058Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in 25mM Tris-Hcl pH 7.2, 1mM EDTA and 50% glycerol. Reconstitute the protein with 100 µl of Millipore water.
Protein	The E.Coli derived recombinant protein contains the CMV Pp28 (UL99) immunodominant regions, 130-160 amino acids.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.