


NS3 (HCV/Genotype-1a/Con)



Product Name	NS3 (HCV/Genotype-1a/Con)
Catalog No	IT-004-066Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS. Reconstitute the protein using 100 µl Ultrapure water.
Protein	GST-fused HCV NS3 Genotype-1a consensus immunodominant regions (aa 1193-1460).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: GST-NS3 (HCV/Genotype-1a/Con)

Protein SEQ:

avdfipvenlettmrspvftdnsspavpqsfvahlhaptgsgkstkvpaaayaagkykvlvlnpsvaatlfgaymskah
gidpnirtgvrvtittgspitystygkfladggcsggaydiiicdehstdatsilgigtvldqaetagarlvvlatatppgsvtvph
nieevalsttgeipfygkaiplevikggrhlfichskkkcdelaaklvalginavayyrgldvsviptsgdvvvvatdalmgtft
gdfdsvidcntcvt