

ORF2 (403-461 a.a.)(HEV/Burma)



Product Name	ORF2 (403-461 a.a.)(HEV/Burma)
Catalog No	IT-004-090Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	N-terminal GST and C-terminal 6xHis tagged the immunodominant regions of ORF2 403-461 a.a. (HEV/Burma) (Genbank# P29326).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Stable for 6 months at -20°C from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

