

ORF3 (82-114 a.a.)(HEV/Burma)



Product Name	ORF3 (82-114 a.a.)(HEV/Burma)
Catalog No	IT-004-093Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	N-terminal GST and C-terminal 6xHis tagged the immunodominant regions of ORF3 82-114 a.a. (HEV/Burma) (Genbank# P69616).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Stable for 6 months at -20°C from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

