



Human TIM-4 (Ectodomain)

Product Name	Human TIM-4 (Ectodomain)
Catalog No	IT-004-TIM-4p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged ectodomain of human TIM-4 protein (a.a. 25-313) (Genebank #: AAH08988). TIM-4 is the member of T-cell immunoglobulin and mucin domain-containing proteins.
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

