



## Nef (HIV-1/Clade A)

---

Product Name	Nef (HIV-1/clade A)
Catalog No	IT-001-013Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS with 8M urea. Reconstitute the protein with 100µl of Millipore water
Protein	C-terminal 6xHis tagged Nef (HIV-1/Clade A)(Genebank accession no. AAC69296.1)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at -80°C; <u>AVOID FREEZING/THAWING</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Nef (HIV-1/clade A) protein

Note: