



## p24 (HIV-1/Clade A)

---

Product Name	p24 (HIV-1/Clade A)
Catalog No	IT-001-016Ep
Source	Recombinant viral protein purified E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS. Add 0.1 ml of ultrapure H <sub>2</sub> O to reconstitute the protein after receiving. Spin briefly before reconstitution.
Protein	C-terminal 6xHis tagged HIV-1/clade A p24 protein (Genebank No. AF004885)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	> 95%
Storage	After reconstitution, store at 4°C for short term. For long term storage, aliquot the protein and store at -80°C. Stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

