




p27(SIV/mac239)

Product Name	p27(SIV/mac239)
Catalog No	IT-001-018p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% NP-40.
Protein	C-terminal 6xHis tagged p27(SIV/mac239) protein (Genebank No. M33262)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: p27(SIV/mac239)

Note: