




p24 (HIV-1/Clade C)

Product Name	p24 (HIV-1/Clade C)
Catalog No	IT-001-019Ep
Source	Recombinant viral protein purified E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS. Add 0.1 ml of ultrapure H ₂ O to reconstitute the protein after receiving. Spin briefly before reconstitution.
Protein	C-terminal 6xHis tagged HIV-1/clade C p24 protein (Genebank No. AY463234)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	> 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E

E: p24 (HIV-1/Clade C)

Note: