




gp120(MG505.W0M.ENV.A2)

Product Name	gp120(MG505.W0M.ENV.A2)(Clade A)
Catalog No	IT-001-101p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120 (MG505.W0M.ENV.A2)(Clade A) protein (a.a.30-506) (Genebank No.DQ208449).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>90%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: (MG505.W0M.ENV.A2)(Clade A) protein

Note: