



gp120(p6244_13.B5.4576)(HIV-1)

Product Name	gp120(p6244_13.B5.4576)(Clade B)
Catalog No	IT-001-127p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120 (p6244_13.B5.4576)(Clade B) protein (a.a.30-523) (Genebank No.EU289191).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(p6244_13.B5.4576) protein

Note: