



gp120(p9014_01.TB1.4769)(HIV-1)

Product Name	gp120(p9014_01.TB1.4769)(Clade B)
Catalog No	IT-001-131p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120 (p9014_01.TB1.4769)(Clade B) protein (a.a.31-530) (Genebank No.EU289195).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(p9014_01.TB1.4769) protein

Note: