



# gp120(p9021\_14.B2.4571)(HIV-1)

---

Product Name	gp120(p9021_14.B2.4571)(Clade B)
Catalog No	IT-001-132p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120 (p9021_14.B2.4571)(Clade B) protein (a.a.32-516) (Genebank No.EU289196).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(p9021\_14.B2.4571) protein

Note: