

gp120(pPRB931_06.TC3.4930) (HIV-1)



Product Name	gp120(pPRB931_06.TC3.4930)(Clade B)
Catalog No	IT-001-134p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120 (pPRB931_06.TC3.4930)(Clade B) protein (a.a.31-497) (Genebank No.EU289198).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(pPRB931_06.TC3.4930) protein

Note: