



Biotinylated gp120(SIVsmE660)

Product Name	Biotinylated gp120(SIVsmE660)
Catalog No	IT-001-142p-Biotin
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of biotinylated protein (0.5 mg/ml) in PBS+40% glycerol
Protein	C-terminal 6xHis tagged gp120 (SIVsmE660) protein (a.a.18-533) (Genebank #: ACQ89608)
Applications	antibody ELISA, protein labeling, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

