



# gp120(SIVcpz MB66)

---

Product Name	gp120(SIVcpz MB66)
Catalog No	IT-001-158p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged gp120(SIVcpz MB66) protein (a.a.39-524) (Genebank No.Q1A261).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(SIVcpz MB66)

Note: