

# gp140(TV-1)(HIV-1/Clade C)



---

Product Name	gp140(TV-1)(HIV-1/Clade C)
Catalog No	IT-001-169p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged gp140(TV-1)(HIV-1/Clade C) (a.a.30-687) (Genebank #: ACF75938)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp140(TV-1)(HIV-1/Clade C)