

# gp120(Q461.d1)(Clade A)



---

Product Name	gp120(Q461.d1)(Clade A)
Catalog No	IT-001-175p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged gp120(Q461.d1)(Clade A)(a.a.34-512) (Genebank Accession No. AF407155_1)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

