

gp120(BG505.W6M.ENV.C1.A204E) (Clade A)



Product Name	gp120(BG505.W6M.ENV.C1.A204E)(Clade A)
Catalog No	IT-001-177p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120(BG505.W6M.ENV.C1)(Clade A)(IT-001-176p) A204E mutant (a.a.30-508) (Genbank No. ABA61515)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(BG505.W6M.ENV.C1)(Clade A)
A204E protein

Note: