




gp120(Con_B)(HIV-1)

| | |
|----------------|---|
| Product Name | gp120(Con_B) |
| Catalog No | IT-001-CONBp |
| Source | Viral protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS. |
| Protein | 6xHis tagged HIV-1 gp120(Con_B/2006) protein (a.a.34-518) (http://www.hiv.lanl.gov/cgi-bin/CONSENSUS_DOWNLOAD/ConsensusDownloader.cgi). |
| Applications | WB standard, antibody ELISA, etc |
| Purity | >= 95% |
| Storage | Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: gp120(Con_B) protein

Note: