



# gp120(AC10.0)(HIV-1)

---

|                |  |
|----------------|--|
| Product Name   | gp120(AC10.0)  |
| Catalog No     | IT-001-RB5p  |
| Source         | Viral protein purified from 293 cell culture   |
| Formulation    | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.   |
| Protein        | 6xHis tagged HIV-1 gp120 (AC10.0)(Clade B) protein (a.a.36-514) (Genebank No. AY835446).   |
| Applications   | WB standard, antibody ELISA, etc   |
| Purity         | ≥ 95%  |
| Storage        | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |

