




gp120(pTRJO4551)(HIV-1)

| | |
|----------------|--|
| Product Name | gp120(pTRJO4551)(Clade B) |
| Catalog No | IT-001-RB7p |
| Source | Viral protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing containing 0.1% BSA and 25% glycerol. |
| Protein | 6xHis tagged HIV-1 gp120 (pTRJO4551)(Clade B) protein (a.a.38-512) (Genebank No. AY835450). |
| Applications | WB standard, antibody ELISA, etc |
| Purity | >95% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: gp120(pTRJO4551) protein

Note: