




gp120(Du156.12)(HIV-1)

Product Name	gp120(Du156.12)(Clade C)
Catalog No	IT-001-RC1p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 25% glycerol.
Protein	6xHis tagged HIV-1 gp120 (Du156.12)(Clade C) protein (a.a.40-505) (Genebank No. DQ411582).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(Du156.12) protein

Note: