




# gp120(ZM233M.PB6)(HIV-1)

---

Product Name	gp120(ZM233M.PB6) (Clade C)
Catalog No	IT-001-RC6p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 25% glycerol.
Protein	6xHis tagged HIV-1 gp120 (ZM233M.PB6)(Clade C) protein (a.a.41-490) (Genebank No. DQ388517).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(ZM233M.PB6) protein

Note: