


HA1(A/Christchurch/16/2010) (H1N1)



Product Name	HA1(A/Christchurch/16/2010)(H1N1)
Catalog No	IT-003-00103p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA(A/Christchurch/16/2010)(H1N1) (a.a.18-344) (See a.a. sequence below)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(A/Christchurch/16/2010)(H1N1)

Protein Sequence:

DTLCIGYHANNSTDTVDTVLEKNVTVTHSVNLLEDKHNGKLCCKLRGVAPLHLGKCNIAGWILGN
PECESLSTASSWSYIVETSSSDNGTCYPGNFIDYEELREQLSSVSSFERFEIFPKTSSWPDHDSNKGV
TAACPHAGAKSFYKNLIWLVKKGNSYPTLSKSYINDKGKEVLVLWGIHHPSTSADQQSLYQNAD
AYVFGTSRYSKFKPEIAIRPKVRNQEGRMNYYWTLVEPGDKITFEATGNLVAPRYAFAMERNA
GSGHISDTPVHDCNTTCQTPKGAINSLPFQNIHPITIGKCPKYVKSTKLRLATGLRNVPSIQSR