

HA(Head Region)(A/Brisbane/59/2007)



Product Name	HA(Head region)(A/Brisbane/59/2007)
Catalog No	IT-003-0012Hp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Head region)(A/Brisbane/59/2007) protein (a.a.51-301) (Genebank #: ACA28846)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

250 kD
150 kD
100 kD
75 kD
50 kD
37 kD
25 kD
20 kD
15 kD
10 kD

E: HA(Head region)(A/Brisbane/59/2007)