

HA₁ (A/Memphis/1/1987(H₁N₁))



Product Name	HA ₁ (A/Memphis/1/1987(H ₁ N ₁))
Catalog No	IT-003-0019p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA ₁ (A/Memphis/1/1987(H ₁ N ₁)) protein (a.a.18-345) (Genebank #: ABQ44416)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(A/Memphis/1/1987(H₁N₁))